


Additional Biosignal Sensor Information

Foot Switch Sensor



DS-X01 Bagnoli System Package

SP-X06	Pre-Amplifier for Bagnoli System
DC-F01	0.5" round FSR

Specifications

Range	Resolution	Bandwidth
0-100 p.s.i.	1.0%	0-500 Hz

DS-X02 Myomonitor System Package

SP-X13	Pre-Amplifier for Myomonitor System
DC-F01	0.5" round FSR

Optional Components

DC-F02	0.2" round FSR
DC-F03	1.5" square FSR

Goniometer Sensor



DS-X03 Bagnoli System Package

SP-X02	Goniometer for Bagnoli System
SC-K02	Goniometer Mounting Kit

Specifications

Range	Resolution	Bandwidth
± 180°	1.0°	0-500 Hz

DS-X04 Myomonitor System Package

SP-X09	Goniometer for Myomonitor System
SC-K02	Goniometer Mounting Kit

EKG Sensor



DS-X09 Bagnoli System Package

SP-X07	EKG Sensor for Bagnoli System
SC-K04	EKG Sensor Kit


Specifications

Range	Resolution	Bandwidth
5 mV _(Pk-Pk)	4 uV _(Pk-Pk)	0.5-30 Hz

DS-X10 Myomonitor System Package

SP-X14	EKG Sensor for Myomonitor System
SC-K04	EKG Sensor Kit

Respiration Sensor



DS-X11 Bagnoli System Package

SP-X17	Pre-Amplifier for Bagnoli System
DC-F04	Respiration Sensor Strap

Specifications

Range	Resolution	Bandwidth
0- 3/4 in.	1/16 in.	0-2 Hz

DS-X12 Myomonitor System Package

SP-X18	Pre-Amplifier for Myomonitor System
DC-F04	Respiration Sensor Strap

Multi-Axis Accelerometer Sensor



DS-X07 Bagnoli System Package

SP-X04	2-Axis Accelerometer for Bagnoli System

Specifications

Range	Resolution	Bandwidth
± 2 g	0.037 g	0-2000 Hz

DS-X08 Myomonitor System Package

SP-X11	2-Axis Accelerometer for Myomonitor System